Morningside Ministries

Skill Checklists to Accompany Taylor's Clinical Nursing Skills: A Nursing Process Approach, 2nd Edition Wolters Kluwer/Lippincott Williams & Wilkins

Name: _____

Date: _____

Unit: _____

Position:

Instructor/Evaluator:

Position:

		Testing Stool for Occult Blood	
Met	Unmet	Goal: An uncontaminated stool sample is obtained following collection guidelines and transported to the laboratory within the recommended time frame, without adverse effect.	Comments
		1. Identify the resident. Discuss with resident the need for a stool sample. Explain to resident the process by which the stool will be collected, either from a bedpan, commode, or plastic receptacle in toilet.	
		2. Close curtains around bed or close door to room if possible.	
		3. Perform hand hygiene.	
		4. Place the plastic collection receptacle in the toilet, if applicable. Assist the resident to the bathroom or onto the bedside commode, or assist the resident onto the bedpan. Instruct resident not to urinate or discard toilet paper with the stool, which may contaminate the specimen.	
		5. After the resident defecates, assist the resident out of the bathroom, off the commode, or remove the bedpan. Perform hand hygiene and put on disposable gloves.	
		6. With wooden applicator, apply a small amount of stool from the center of the bowel movement onto one window of Hemoccult testing card. With opposite end of wooden applicator, obtain another sample of stool from another area and apply a small amount of stool onto second window of Hemoccult card.	
		7. Close flap over stool samples.	
		8. If testing at bedside, open flap on opposite side of card and place two drops of developer over each window and wait the time stated in the manufacturer's instructions.	
		9. Observe card for any blue areas.	
		10. Discard Hemoccult testing slide appropriately, according to facility policy. Remove gloves and perform hand hygiene.	